

Crestpoint ROUNDS™

NOVEMBER
2019
Volume 3
Issue 3

PUBLISHED FOR ALL OPHTHALMIC AUDIENCES: SURGEONS, YOUNG EYE SURGEONS, NURSES/TECHNICIANS AND PURCHASING TEAMS

Special Edition ~ In Memory of Dr. Spencer Thornton

CRESTPOINT—D&K—GEUDER—MANI—MR.PHACO/MR.FEMTO—STERIMEDIX—QINSTRUMENTS—AMCON



IN THIS ISSUE

Duckworth & Kent	Cover, 2, 3
Mr. Phaco/Mr. Femto	4

In Memory of Dr. Spencer Thornton


Over the past 40 years, Dr. Thornton – a true innovator, designed and developed product with Crestpoint, Duckworth & Kent and the entire ophthalmic industry.

This Edition of ROUNDS is dedicated to our innovator, teacher, mentor, author, magician, actor, and our dearest friend. With this – please carry on innovation and creativity to continue the legacy of our dear colleague.

Sincerely,





[Click here for our new Catalog!](#)



9-510R Fine Thornton Fixation Ring

3/4 ring, 13.0mm, 9 point fixation, pivot swivel, round handle.






Dr. S. Thornton's finest innovations


Duckworth & Kent are proud to have had the opportunity to work with Dr. S. Thornton over many years to design and produce an exquisite array of reusable ophthalmic instrumentation.

Thornton Fixation Rings


A varied range of precision instruments to suit all surgeons needs
Flat handle, round handle, angled, pivot swivel.




9-515R



9-510 & 9-510R




9-509




Thornton Radial Markers


First-rate instruments for accurate markings
Centre pointed, round handle.




9-730




9-732



9-853




9-733




9-734

Thornton Diamond Knives


3 distinctive blade variations
Single edge, triple edge, micrometer handle.




5-330-1



5-362




5-360-1



Thornton Comfort Speculum

Finer, delicate speculum blades provide comfort for patient
Functionally designed with open blades to minimize interference with
intraocular instrumentation.



9-581F




D&K Capsulorhexis Forceps

AAO 2019 Most Popular

2-716G-10RE Jones - Inamura Round Handle Capsulorhexis Forceps

Now available to view in 360°
Visit - www.duckworth-and-kent.com




E-Range



2-716G-9E Calladine - Inamura Flat Handle Capsulorhexis Forceps

Now available to view in 360°
Visit - www.duckworth-and-kent.com




E-Range




D&K Instruments

AAO 2019 Most Popular


1-110 Vannas Scissors (straight)




1-705 Osher IOL Cutter (angled)




2-500E McPherson Tying Forceps



E-Range




2-110E Bonn Suturing Forceps (0.12mm)




E-Range



8-657 DK Aspiration Handpiece



6-250 DK Sinsky Hook



E-Range Instruments


Elegant - Ergonomic - Precision Engineered

This range represents the application of cutting edge technology fused with over 60 years of engineering knowledge. They offer surgeons an elegant, well balanced and ergonomically designed handle, providing our consistent high quality.

Design registration number 004383396-0001



3-205 DK Needle Holder (curved)



AAO 2019

Thank You

from all of us at
Duckworth & Kent and Crestpoint
We have had another successful
AAO!

We would like to thank all those who
visited us at our booth, and purchased
our reusable titanium instruments

We also hope you found time to enjoy the
many sights of San Francisco.





Mr.Phaco® / Mr.Femto®



One Source service and support for your ophthalmic needs. New innovative instruments, repair/replacement program, custom instruments. Excellent reputation for service and quality.

Arun C. Gulani, MD

PH6013 Gulani Harpoon Universal Chopper Manipulator

Double Ended, designed to assist surgeon with manipulation/rotation of cataract and nucleus. Harpoon end 60° curved, tip 0.75mm x 0.51mm x 2.0mm. Chopper end 45° angled shaft, 90° angled tip.

Brochure available!



[Click here for a demonstration with Dr. Gulani!](#)



Arun C. Gulani, MD & Spencer P. Thornton, MD



Henry M. Mitchell, J. Michelle Glossip & Spencer P. Thornton, MD





Robert H. Osher, MD & Spencer P. Thornton, MD

Robert H. Osher, MD

PH9304 Osher Capsulorhexis Ring Marker

Single Ended, 5.75mm inner diameter, 6.44mm outer diameter, round handle. Designed to leave a faint circular imprint on the corneal epithelium to assist the surgeon with precise sizing of the capsulorhexis.

Brochure available!



Contact Us



Our top priority is superb customer service to you.



Crestpoint Management
Crestpoint Ophthalmics
 9100 Watson Road
 St. Louis, MO 63126 USA

GSA/FSS Contract
 V797P-4438B
 Ophthalmic Instruments
 SIN A8a & A8b
 Effective
 5/1/2011 to 4/30/2021

Published by Crestpoint Management. Manufacturer trademarks, part data, and copyrighted materials are property of their respective owners and are used with permission. D&K, Geuder, MANI product ranges in this publication are exclusive in the USA / NA.

- Phone: (800) 756-7773
- Phone: (314) 849-7773
- Fax: (314) 849-7776
- E-mail: custsvc@crestpointmgt.com
- Web: www.crestpointmgt.com

Copyright 2019 Crestpoint Management / Crestpoint Ophthalmics
 All rights reserved. All information in this publication is subject to change without notice.

ZCC483NL (11/2019)