

# Crestpoint ROUNDS™

JANUARY  
2018  
Volume 2  
Issue 1

PUBLISHED FOR ALL OPHTHALMIC AUDIENCES: SURGEONS, YOUNG EYE SURGEONS, NURSES/TECHNICIANS AND PURCHASING TEAMS

## IN THIS ISSUE

### IRRIGATION AND ASPIRATION

GEUDER.....	1, 3, 4
Duckworth & Kent.....	2
MANI.....	2
STERIMEDIX.....	3
Mr.Phaco/Mr.Femto.....	4
GEUDER QUALITY.....	4
Crestpoint Contact Information.....	4

One  
Source®

## Irrigation/Aspiration Handpieces (I/A)

We appreciate the opportunity to serve you:

This issue of *Crestpoint ROUNDS* is dedicated to Irrigation/Aspiration Handpieces (I/A) – **featuring GEUDER.**

**Crestpoint Management is the exclusive USA Distributor for GEUDER Handheld Instruments.**

GEUDER AG is one of the leading manufacturers of ophthalmic surgical instruments and systems in Germany and abroad. With extensive involvement in research and development, the GEUDER Group sets standards and demonstrates a future-oriented company philosophy.

Exporting to more than 90 countries, GEUDER supplies a complete range of over 3,000 precision surgical instruments for anterior and posterior eye segment surgery.

Look inside to see how we focus on **excellent product quality and customer service – dedicated to you, and your organization.**

Sincerely,  
Your Crestpoint Team


J. Michelle Glossip  
Chief Executive

Henry M. Mitchell  
President

Dot C. Smith  
Customer Resource  
Director

Michael Dill  
Customer Resource  
Manager





# Duckworth & Kent


Original manufacturer of reusable titanium surgical instruments for ophthalmology.  
[Click here to download the D&K Irrigation/Aspiration PDF Catalog](#)

- 8-711N - Single Thread
- 8-711NL - Single Thread with Luer Lock
- 8-711-1N - Multi Thread
- 8-711-1NL - Multi Thread with Luer Lock



## I/A

- Round handle, length 115mm
- Used with DK I/A tips
- I/A tips purchased separately
- Single/Multi thread/Luer Lock



- 23 gauge, 0.65mm tube diameter
- Two 0.4mm irrigation ports
- Curved shaft, tube length 15mm
- Round long handle, length 107mm - 8-652
- Round short handle, length 80mm - 8-652S



- 23 gauge, 0.65mm tube diameter with textured tip
- 0.3mm aspiration port
- Curved shaft, tube length 11mm
- Round long handle, length 103mm - 8-657
- Round short handle, length 77mm - 8-657S



# MANI®



Single use Knives, Sutures, Trocars for all ophthalmic procedures. The crosslapped edge of MANI®'s original hard fiber stainless steel produces the "Ultimate Sharpness."

### Short-Head Slit Knife

Main Application: Corneal incision, Scleral incision, 6pcs/box, disposable



Effective for a patient who has sunken / deep-set eye (enophthalmos) or narrow lid fissure.



Position of knife angle moved 2mm towards the tip and the blade became smaller with extreme sharpness as our conventional models.

### Product Name

- Slit Angled 2.0mm (Short-Head)
- Slit Angled 2.2mm (Short-Head)
- Slit Angled 2.4mm (Short-Head)
- Slit Angled 2.75mm (Short-Head)
- Slit Angled 3.0mm (Short-Head)

### Order #

- MSL20SH
- MSL22SH
- MSL24SH
- MSL27SH
- MSL30SH




# GEUDER®

Leading manufacturer of high quality surgical instruments for ophthalmology. It has always been the maxim of GEUDER to guarantee the highest degree of precision and reliability.


### D. D. KOCH - T. KOHNEN Aspiration Handpiece

- G-22100**  
**G-22102**  
0.3 mm aspiration port  
textured tip
- G-22100** 22 gauge / 0.7 mm
- G-22102** 22 gauge / 0.6 mm




### D. D. KOCH - T. KOHNEN Irrigation Handpiece

- G-22101**  
front opening  
22 gauge / 0.7 mm



### Aspiration Handpiece

- G-22105**  
**G-22106**  
0.25 mm aspiration port  
conical tip  
oval tube 0.95 x 0.55 mm
- G-22105** smooth
- G-22106** textured




### SCHMACK I/A Handpiece

- G-22632**  
**G-22633**  
21 gauge / 0.8 mm
- G-22632** irrigation  
0.5 mm port
- G-22633** aspiration  
0.3 mm port




### Aspiration Handpiece

- G-32011**  
0.3 mm aspiration port  
textured, conical tip  
21 gauge / 0.8 mm



### Irrigation Handpiece

- G-32012**  
dual 0.5 mm irrigation ports  
textured, conical tip  
21 gauge / 0.8 mm




# Sterimedix®



Complete range of single use devices for all ophthalmic surgery procedures.


### Irrigating/Aspirating Cannulas and Handpieces

Code	Gauge	Description
<b>M5570</b>	24G	Soft Outer I/A Handpiece for 2.2mm incisions. Straight tip.
<b>M5573</b>	24G	Soft Outer I/A Handpiece for 2.2mm incisions. 90° angled tip.



### Irrigating/Aspirating Cannulas and Handpieces

Code	Gauge and size	Description
<b>M5531-set</b>	21G x 5/8" 0.81 x 16mm	Bimanual Irrigation Handpiece. Thin Wall. Curved with twin 0.5mm side ports.
	21G x 5/8" 0.81 x 16mm	Bimanual Aspiration Handpiece. Curved with 0.4mm side port and micro-etched tip.
<b>M5331-set</b>	21G x 5/8" 0.81 x 16mm	Bimanual Irrigation Handpiece. Curved with 2 x 0.4mm irrigation ports and tapered tip.
	21G x 5/8" 0.81 x 16mm	Bimanual Aspiration Handpiece. Curved with 0.3mm aspiration port and micro-etched tapered tip.
<b>M5535-set</b>	23G x 5/8" 0.64 x 16mm	Bimanual Irrigation Handpiece. Thin Wall. Curved with twin 0.4mm side ports.
	23G x 5/8" 0.64 x 16mm	Bimanual Aspiration Handpiece. Curved with 0.3mm side port and micro-etched tip.





# Mr.Phaco® / Mr.Femto®

One Source service and support for your ophthalmic needs. New innovative instruments, repair/replacement program, custom instruments. Excellent reputation for service and quality.

### I/A Micro Tip Cannula

**MMP250**  
21 gauge, 0.25mm,  
straight shaft,  
stainless steel, reusable



--  
\*20 gauge, 0.50mm,  
straight shaft,  
specify **MMP257**

### I/A Micro Tip Cannula

**MMP252**  
21 gauge, 0.25mm,  
angled 45°,  
stainless steel, reusable



--  
\*21 gauge, 0.30mm,  
angled 45°,  
specify **MMP260**

### I/A Micro Tip Cannula

**MMP262**  
21 gauge, 0.3mm,  
curved tip,  
stainless steel, reusable



### I/A Micro Tip Cannula

**MMP251**  
21 gauge, 0.25mm,  
angled 25°,  
stainless steel, reusable



### I/A Micro Tip Cannula

**MMP253**  
21 gauge, 0.25mm,  
angled 90°,  
stainless steel, reusable



--  
\*21 gauge, 0.30mm,  
angled 90°,  
specify **MMP261**

### I/A Handpiece


**PH2003**  
(tips purchased separately)  
stainless steel, reusable



## GEUDER QUALITY



The company's principles are to guarantee maximum precision and reliability, combining technological achievements with a tradition of skilled craftsmanship.



Based in Heidelberg, Germany, the company supplies a complete range of surgical instruments used for anterior and posterior eye segment surgery. Through intensive contact and cooperation with experienced surgeons specialized in ophthalmic procedures, GEUDER

is able to develop continuously innovative eye surgery products which surpass international standards and are highly appreciated. Due to extensive involvement in research and development, the GEUDER Group sets standards with a future-oriented company philosophy.

## ABOUT US

### Mission Statement

Crestpoint Ophthalmics is dedicated to being the foremost, industry leading, single-source provider of the world's highest quality ophthalmic products.

We're your source for 5,000+ hand-held instruments, instrument repair service, repair/replacement program, custom instruments and reusable/disposable instruments. Our top priority is superb customer service to you.



All information in this publication is subject to change without notice.

Copyright 2018 Crestpoint Ophthalmics. All rights reserved.

Trademarks and copyrighted materials are property of their respective owners and are used with permission.



All product ranges in this Newsletter are exclusive in the USA.

Duckworth & Kent products are exclusive in the USA and Canada.

## CONTACT US

### Crestpoint Management

9100 Watson Road  
St. Louis, MO 63126 USA

☎ Phone: (800) 756-7773

☎ Phone: (314) 849-7773

☎ Fax: (314) 849-7776

✉ E-mail: [custsvc@crestpointmgt.com](mailto:custsvc@crestpointmgt.com)

🌐 Web: [www.crestpointmgt.com](http://www.crestpointmgt.com)

ZCC483NL (01/2018)