

Crestpoint ROUNDS™

APRIL
2018
Volume 2
Issue 3

PUBLISHED FOR ALL OPHTHALMIC AUDIENCES: SURGEONS, YOUNG EYE SURGEONS, NURSES/TECHNICIANS AND PURCHASING TEAMS

IN THIS ISSUE

Duckworth & Kent.....	Cover, 2, 3
GEUDER.....	Cover, 3, 4
Mr.Phaco./Mr.Femto.....	4

Looking forward to seeing you at the
ASCRS-ASOA 2018 Meeting
Visit Us At Booths 2227 & 2324

**One
Source**
SERVICE AND SUPPORT

Forceps

Greetings Ophthalmic Community:

In this issue we are highlighting **forceps** as the first in a series of newsletters that are focused on specific instrument categories — featuring the finest quality surgical instruments from **Duckworth & Kent®** and **Geuder®**.

Crestpoint Ophthalmics is the exclusive USA Distributor for

- **Duckworth & Kent® Titanium Surgical Instruments**
- **GEUDER® Handheld Instruments**
- **Mr.Phaco® / Mr.Femto®**
- **MANI® Knives and Sutures**
- **Sterimedix®**

Crestpoint ROUNDS is published monthly for all ophthalmic audiences: surgeons, young eye surgeons, nurses/technicians, and purchasing teams — and our mission is to keep the Ophthalmic Community informed by highlighting the **best of the best** product solutions that enhance and improve the efficiency of ophthalmic surgical techniques.



Duckworth & Kent
titanium surgical instrument manufacturer

Innovation is at the heart of our philosophy and is driven by our enduring commitment to provide surgeons with only the highest quality reusable titanium ophthalmic instrumentation.

We are proud to announce and unveil our latest line of instrumentation, the "E-Range".

This range represents the application of cutting edge technology fused with over 50 years of engineering knowledge. These instruments offer surgeons an elegant, well balanced and ergonomically designed handle whilst still providing our high quality, consistent and durable functionality.



E-RANGE
Instruments


Geuder®

Precision made in Germany

GEUDER AG is one of the leading manufacturers of ophthalmic surgical instruments and systems in Germany and abroad. Through their extensive work in research and development, the GEUDER Group sets standards and demonstrates a future-oriented company philosophy.

Exporting to more than 90 countries, GEUDER supplies a complete range of over 3,000 precision surgical instruments for anterior and posterior eye segment surgery.






Duckworth & Kent

Original manufacturer of titanium surgical instruments for ophthalmology.

Capsulorhexis Forceps




- Sharp pointed serrated tips
- Cross action tips, 1.5mm width at pivot box
- Curved shaft, tip to pivot point 8.5mm / 9.5mm / 10.5mm
- Marks on shaft 2.5mm and 5mm denote desired size of capsulorhexis
- Tips angled 45° from shaft



Product No.	Product Name	Variation
2-716GER8	Inamura Round Handle Capsulorhexis Forceps	Pivot Point 9.5mm, long round handle 110mm
2-716GE-1	Inamura Flat Handle Capsulorhexis Forceps	Pivot Point 10.5mm, long flat handle 111mm
2-716GE-1S	Inamura Flat Handle Capsulorhexis Forceps	Pivot Point 10.5mm, short flat handle 86mm
2-716GER8-1	Inamura Round Handle Capsulorhexis Forceps	Pivot Point 10.5mm, long round handle 111mm
2-716GER8-1S	Inamura Round Handle Capsulorhexis Forceps	Pivot Point 10.5mm, short round handle 86mm
2-716GER8-2	Inamura Round Handle Capsulorhexis Forceps	Pivot Point 8.5mm, long round handle 110mm


Long Handle Forceps

Handle length 115mm




2-104E

Notched




2-510-1E

Utility



2-114E
2-116E
2-135E


Toothed



2-504E
2-504ER8 Round handle
2-527E


Tying

Short Handle Forceps




2-100E

Notched




2-500E
2-522E

Tying




2-110-3E

Toothed




2-500-4E

Utility



2-110E
2-118E
2-132E


Toothed



2-2195E

Cilia

Handle length 90mm



E-RANGE
Instruments




GEUDER®

Leading manufacturer of high quality surgical instruments for ophthalmology. It has always been the maxim of GEUDER to guarantee the highest degree of precision and reliability.

Capsulorhexis Forceps

G31299



Angled, 12.5mm platforms, cross serrated tips

G31305



Curved, serrated tips

G32944

HATTENBACH



Angled, mini platforms, 25gauge/0.5mm

G31304



Curved, serrated delicate tips – specify **G31304**
Angled, serrated heavy tips – specify **G31307**

G32928

INAMURA



Cross-action, serrated tips

G32926

AUFFARTH



Flat handle, cross-action, for incision size 1.6mm, total length 100mm

Colibri Forceps

G18950



Tying platform, inter-meshing teeth – specify **G18950**
Micro, with tying platform – specify **G18960**

Dressing Forceps

G18792



1x2 meeting teeth, 1.0mm

Pre-Chop Forceps

G31519

AKAHOSHI



Cross-action

G31724

AKAHOSHI



Micro-incisions, cross-action, straight, round handle, length 120mm

Suturing / Tying Forceps

G19001

CASTROVIEJO



Straight, 0.28mm, with tying platform

G19022

CASTROVIEJO



Straight, with tying platform, 0.25mm – specify **G19022**
Straight, with tying platform, 0.3mm – specify **G19023**
Straight, with tying platform, 0.5mm – specify **G19024**

Tying / Fragment Forceps

G19030

TUEBINGEN



Straight – specify **G19030**
Angled – specify **G19032**
Curved – specify **G19040**

G31624

GASKIN



Angled, 11.0mm

G31623

KELMAN-MCPHERSON



Angled, 7.5mm

G32947

HATTENBACH



Angled, straight platform, 25 gauge/0.5mm



Mr.Phaco® / Mr.Femto®

One Source service and support for your ophthalmic needs. New innovative instruments, repair/replacement program, custom instruments. Excellent reputation for service and quality.

Single Use Forceps

Q instruments

Innovative manufacturing technology.
Precise, consistent, ergonomic.

- Efficient utilization and safety
- Logistical cost efficiency
- Guaranteed performance and sterility



G10010 Capsulorhexis Forceps

Sharp tips, 11.5mm curved jaws
80mm handle, overall length 102mm



G10020 Suturing Forceps

0.12mm 1x2 teeth, 6.5mm tying platforms
80mm handle, overall length 100mm



G10031 Tying Forceps, angled

6.5mm tying platforms
80mm handle, overall length 100mm



G10032 Notched Forceps

10mm tips, 0.25mm notch
80mm handle, overall length 100mm



G10050 IOL Forceps, angled

0.5mm x 8.7mm tips, 45° angled shaft
80mm handle, overall length 100mm



duckworth & kent
titanium surgical instrument manufacturer

We have a wide range of precision engineered forceps.

For more information please visit www.duckworth-and-kent.com

SEE YOU IN WASHINGTON D.C.!

Looking forward to seeing you in April at the ASCRS-ASOA 2018 American Society of Cataract and Refractive Surgery Meeting

April 13-17, 2018

Walter E. Washington Convention Center

Visit Us At Booths 2227 & 2324

Geuder® Precision made in Germany

Crestpoint MANAGEMENT
...superior products...no boundaries...

Spring / Summer 2018 FEATURED INSTRUMENTS

This catalog is designed to facilitate your choices among our 5,000+ product range.

Precision, Consistency, Excellence. Finest Materials Used in Manufacturing for YOUR Procedures. SEE for Yourself!

We Proudly Supply the Finest and Very Best.

Thank you for Appreciating OUR Service and the Quality of OUR Product Range.

Crestpoint Management
Exclusive USA Distributor

ABOUT US

Mission Statement

Crestpoint Ophthalmics is dedicated to being the foremost, industry leading, single-source provider of the world's highest quality ophthalmic products.

We're your source for 5,000+ hand-held instruments, instrument repair service, repair/replacement program, custom instruments and reusable/disposable instruments. Our top priority is superb customer service to you.



GSA/FSS Contract
V797P-4438B
Ophthalmic Instruments
Surgical Devices
SIN A8a & A8b
Effective
5/1/2011 to 4/30/2021

All information in this publication is subject to change without notice. Trademarks and copyrighted materials are property of their respective owners and are used with permission.

NWBOC
WVOSB CERTIFIED-WOMAN OWNED SMALL BUSINESS

All product ranges in this Newsletter are exclusive in the USA.

Duckworth & Kent products are exclusive in the USA and Canada.

Copyright 2018 Crestpoint Ophthalmics. All rights reserved.

CONTACT US

Crestpoint Ophthalmics
9100 Watson Road
St. Louis, MO 63126 USA

Phone: (800) 756-7773

Phone: (314) 849-7773

Fax: (314) 849-7776

E-mail: custsvc@crestpointmgt.com

Web: www.crestpointmgt.com

ZCC483NL (04/2018)