

# Crestpoint ROUNDS™

APRIL  
2019  
Volume 3  
Issue 2

PUBLISHED FOR ALL OPHTHALMIC AUDIENCES: SURGEONS, YOUNG EYE SURGEONS, NURSES/TECHNICIANS AND PURCHASING TEAMS

## Special Edition ~ See You In San Diego!

WELCOME TO CRESTPOINT—D&K—GEUDER—MANI—MR.PHACO/MR.FEMTO—STERIMEDIX—QINSTRUMENTS—AMCON




### DUCKWORTH & KENT

**Page 2**  
World leader in reusable titanium surgical instruments for ophthalmology.



### GEUDER

**Page 3**  
High quality stainless steel surgical instruments for ophthalmic surgery.



### MANI

**Page 3**  
MANI's fine processing technique creates the crosslapped edge (patented).



### MR.PHACO

**Page 4**  
Innovative and custom instruments, repair/replacement program.



### STERIMEDIX

**Page 2**  
Single-use disposable products deliver innovation in ophthalmic surgery.



### Q INSTRUMENTS

**Page 4**  
Single-use ophthalmic instruments.



## Anterior Vitrectomy and Complicated Cases

See you in San Diego at ASCRS 2019!

The owner, exclusive supply partners, new representatives, and our specialized customer resources team will be on hand. Dedicated to the ophthalmic community, with **best-of-the-best products and service**. Always count on Crestpoint for trusted, personalized attention.

**OneSource®** for 5,000+ hand-held surgical instruments, sales, service, repair/replacement programs, research & development, design, engineering, manufacturing, and expert technical support. **Our top priority is superb customer service to you.**

Please come and visit us at ASCRS 2019 booth #s 1907 and 2006. Innovation, service, and best prices for you.

Crestpoint ROUNDS — Always featuring best of the best. Keeping you informed, providing product solutions to enhance and improve efficiency of surgical techniques.

See inside functionality and form, precision at its' best.

### Duckworth & Kent Reusable Titanium

8-730 1/A Tips for Use With 8-711N/8-711-IN Handpieces



8-635-3 Ogawa 1/A Cannula



### Geuder Chromium Plated Surface

G19022 *Castronigo* Suturing Forceps



G31305 Capsulorhexis Forceps



### MANI Ultimate Sharpness

MANIMTR23S Trocar Kit 23G.



MANIMST15MK Straight Knife 15 Degree Marking Sideport



### Sterimedix Continuous Development

M4620 Irrigating Vectis



M2657 *Lewicky* Chamber Maintainer



### Mr.Phaco/Mr.Femto Service and Support

MMP265 *Osher* Sweep for Chemosis



PH8002 *Osher* Vectis Cannula on Handle




### Q instruments Efficient and Disposable

G10020 Suturing Forceps



G10031 Tying Forceps, angled





# Duckworth & Kent Titanium. Reusable.

Original manufacturer of titanium surgical instruments for ophthalmology and known throughout the world for quality and innovation. D&K titanium instruments are reusable, lightweight, non-magnetic, and anodized to reduce glare under the microscope.



## Anterior Vitrectomy Instruments

### Ogawa I/A Cannula

8-635-3 Ogawa I/A Cannula


### 22g Cannula with Handle

8-730-2 DK I/A Tip (J-shaped tip)  
8-730 DK I/A Tip (90° angled tip)  
8-711-1NL I/A Handpiece, Luer Lock



### Sub Tenon Set

9-551 Barraquer Adult Speculum  
1-501B DK Westcott Style Stitch Scissors  
2-510-1E Moorfields Utility Forceps



### Needle Holder

3-203 DK Needle Holder (curved)



### Forceps

2-100E Pierser Notched Forceps  
2-500E McPherson Tying Forceps

### Manipulators / Rotators


6-085 Barrett Duo Nucleus Rotator / Manipulator / Splitter  
6-400 Kuglen Iris Hook and IOL Manipulator  
6-250 DK Sinskey Hook



See the full product range at [www.duckworth-and-kent.com](http://www.duckworth-and-kent.com)



## STERIMEDIX® Single-Use Disposable



Safety. Quality. Flexibility. Cost efficiency. STERIMEDIX single-use, disposable surgical products meet the demands of modern ophthalmic surgical techniques.

### M0278 Sub Tenon Cannula

19 gauge x 1", curved with flattened tip.



### M2273A Rycroft Anterior Chamber Maintainer

27 gauge x 7/8", angled 35° 4.5mm from tip.



### M2280 Transzonular Cannula

27 gauge, angled 33° 2.0mm from tip.




### M2657 Lewicky Chamber Maintainer

20 gauge x 9/64", self retaining, silicone tubing included. 23 gauge, specify M2658.



### M5573 Irrigating/Aspirating Handpiece

24 gauge, soft outer I/A handpiece for 2.2mm incisions. 90° angled tip.




### M5535 Bimanual I/A System

Bimanual irrigation handpiece, 23 gauge x 5/8", thin wall, curved, twin 0.40mm sideports. Bimanual aspiration handpiece, 23 gauge x 5/8", thin wall, curved, 0.3mm sideport, micro etched tip.



### M4620 Irrigating Vectis

Irrigating Vectis. 5mm wide, twin ports.





# GEUDER<sup>®</sup> Chromium Plated Stainless Steel. Reusable.

Leading manufacturer of high quality surgical instruments for ophthalmology, made with chromium plated stainless steel. It has always been the maxim of GEUDER to guarantee the superior level of craftsmanship, precision and reliability.

### G1750 Barraquer Needle Holder

Delicate 10.0mm curved jaws, tag spring, without lock, round handle.



### G31305 Capsulorhexis Forceps

Extra delicate interlocking micro serrated tips, 13.0mm curved shafts, direct action flat handle. Strong with heel stop.



### G19751 Westcott Style Stitch Scissors

Sharp tips, curved standard blades, tag spring, flat handle. Blunt tips, curved standard blades, tag spring, flat handle, specify G19750.



### G19022 Castroviejo Suturing Forceps




0.25mm, 0.12 teeth, tying platforms, straight shafts, flat handle.

### G22100 D.D. Koch - T. Kohnen Aspiration Handpiece



22 gauge/0.7mm, 0.3mm port, textured tip, curved shaft, male connector.

### G32996 Snyder Osher IOL Cutter (for Snyder Osher Forceps G32997)



Long blade 4.5mm, 18 gauge for cutting acrylic and Silicone foldable IOL's. Serrated tips prevent the "watermelon seed" phenomenon, especially vexing with smooth scissors edges and implants. Scissors can section an IOL either in pieces or for hinged folding, facilitating explantation through a small incision. The rounded tips reduce the chances of inadvertently catching iris tissue or capsule, while the fine serrations can help prevent the IOL from slipping out of the jaws during cutting. Strong enough for acrylic IOLs. Short 3.5mm blades, specify G32998.

### G32997 Snyder Osher Forceps (for Snyder Osher IOL Cutter)



(Incision Size 1.9mm) For use with G32996 and G32998. These unique forceps have square-tipped ends with teeth that can firmly grip an IOL or IOL fragment. The three-dimensional gripping prevents the piece from torquing or tilting with the forceps. Small profile instrument fits easily through a paracentesis and a 2.2mm incision. Rounded tips reduce the risk of engaging strands of iris tissue and reduce the risk of damage to the delicate intraocular environment. 2.0mm jaws with 0.5mm serrated tips, 21 gauge lightly curved shaft.




# MANI<sup>®</sup> Single-Use Disposable Ophthalmic Knives

MANI's bevel shape (patented), crosslapped edge, and high-strength stainless steel construction produces the "Ultimate Sharpness" for a true, stable incision with higher accuracy and lower penetration resistance. Color-coded knife handles allow for easy identification. The color is different depending on the tip shape.


### Straight Knife 15 Degree / 22.5 Degree Marking Sideport

MANIMST15MK Straight Knife 15 Degree Marking Sideport. Box of 6  
 MANIMST22MK Straight Knife 22.5 Degree Sideport Marking. Box of 6



### MVR Knife Angled 23G 0.9mm


MANIMVR23A MVR Knife Angled 23G 0.9mm. Box of 6



### Trocar Kit 23G

MANIMTR23S Trocar Kit 23G. Box of 6

Trocar blade in the shape of MVR knife provides ultimate sharpness adopting proven technologies of MANI<sup>®</sup> Ophthalmic Knife. Valved cannula contributes to the simpler operation.



### MVR Knife Series Item Numbers

- |          |        |         |
|----------|--------|---------|
| MVR20ASK | MVR23  | MVR23A  |
| MVR21ASK | MVR24  | MVR24A  |
| MVR23ASK | MVR25  | MVR23AD |
| MVR24ASK | MVR19A | MPK25   |
| MVR19    | MVR20A | MPK27   |
| MVR20    | MVR21A |         |



# Mr. Phaco® / Mr. Femto®

## MMP265 Osher Sweep for Chemosis

Designed to sweep chemosis, subconjunctival anesthesia, BSS, or hemorrhage away from the incision to allow unimpeded wound construction. Single ended, flat handle. Stainless Steel, reusable



VIDEO AVAILABLE

## PH8002 Osher Vectis Cannula on Handle

23 gauge, 3.0mm x 16.0mm, three 0.3mm ports at 11, 12, 1 o'clock, round handle. Cannula only, specify **PH8001**. Miniature Nucleus Vectis is reduced size to allow the surgeon to rescue a nuclear fragment through a small incision. The distal angled ski configuration facilitates nuclear recovery. Malleable tip allows surgeon to quickly configure the optimal shape for nucleus recovery. Stainless Steel, reusable



## H1712 Osher Y Hook

Designed for lens rotation during which the polished dull finger will not mark the lens or tear the capsule if accidental contact occurs. Single ended, straight shaft, ball shaped tip, round handle. 45° angled, specify **H1712A**. Stainless Steel, reusable



# Q instruments® Single-Use Disposable

## G10020 Suturing Forceps

Suturing Forceps  
Stainless Steel tips, disposable



## G10031 Tying Forceps

Tying Forceps, angled  
Stainless Steel tips, disposable



## See You In San Diego at ASCRS Booth #'s 1907-2006, 2007



## Crestpoint Team Spotlight

Stephanie Rena' Gibson



Stephanie is one of our dynamic Customer Resource Team members, working together with our other team members to help ensure that your customer service experience with Crestpoint is always handled in a positive, efficient, and friendly manner.

Stephanie also serves as a contributor to our sales support and expanding marketing operations team, and she is directly involved in working with our product experts to help manage our product catalogs, brochures and literature, and direct mail programs.

We always strive to make your experience the best it can be, and Stephanie is one of our helpful ambassadors. We can't thank our customers enough!

## Contact Us



Our top priority is superb customer service to you.

NWBOC  
WOSB CERTIFIED: WOMAN OWNED SMALL BUSINESS

**Crestpoint Management**  
**Crestpoint Ophthalmics**  
9100 Watson Road  
St. Louis, MO 63126 USA

GSA/FSS Contract  
V797P-4438B  
Ophthalmic Instruments  
Surgical Devices  
SIN A8a & A8b  
Effective  
5/1/2011 to 4/30/2021

Published by Crestpoint Management. Manufacturer trademarks, part data, and copyrighted materials are property of their respective owners and are used with permission. D&K, Geuder, MANI product ranges in this publication are exclusive in the USA / NA.

- Phone: (800) 756-7773
- Phone: (314) 849-7773
- Fax: (314) 849-7776
- E-mail: [custsvc@crestpointmgt.com](mailto:custsvc@crestpointmgt.com)
- Web: [www.crestpointmgt.com](http://www.crestpointmgt.com)

Copyright 2019 Crestpoint Management / Crestpoint Ophthalmics  
All rights reserved. All information in this publication is subject to change without notice.  
ZCC483NL (4/2019)