

Crestpoint ROUNDS™

JUNE
2018
Volume 2
Issue 5

PUBLISHED FOR ALL OPHTHALMIC AUDIENCES: SURGEONS, YOUNG EYE SURGEONS, NURSES/TECHNICIANS AND PURCHASING TEAMS

IN THIS ISSUE

STERIMEDIX®	Cover, 3, 4
Duckworth & Kent®	2
GEUDER®	2

Looking forward to seeing you at
AAO 2018
October 27-30 Chicago

Precision, Single-Use Disposable, and Reusable Surgical Cannulas


Featuring single-use disposable ophthalmic surgical products from STERIMEDIX®

Benefits of Single-Use Devices:

- ✓ Safety
- ✓ Quality
- ✓ Improved patient care
- ✓ Cost efficiency and inventory control
- ✓ Regulatory compliance

User-Friendly Packaging Options:

- ✓ Sterile—packaged product in boxes of 10 or 100 reduces cost and helps ensure a ready supply of safety stock.
- ✓ Non-Sterile—individually labeled with complete traceability for kit packers.



STERIMEDIX® Color-Coded Hubs
Twelve colors provide a user-friendly reference and product key for 16 gauge–30 gauge hubs.

STERIMEDIX® Precision Manufacturing


The requirements of eye surgery place a great responsibility on the manufacturer to maintain the highest standards of quality at all times. The company's quality systems and practices are regularly audited and have earned ISO 13485, ISO 9001, CE and U.S. FDA registration.




Sterile Packaged Product



ZSC369 Catalog




Kit Packers Non-Sterile Bulk Packaging



Duckworth & Kent Titanium-Reusable


Coaxial Handpieces & Tips

Multiple Tip Styles Available



- 22 gauge, 0.7mm tube diameter
- 0.3mm aspiration port
- 90° angled tip

22g 8-731




- 20 gauge, 0.9mm tube diameter
- 0.3mm aspiration port
- 45° angled, textured tip

20g 8-732-7



- 8-711N — Single Thread
- 8-711NL — Single Thread with Luer Lock
- 8-711-1N — Multi Thread
- 8-711-1NL — Multi Thread with Luer Lock



Separates for cleaning internal parts




(Silicone sleeve not included)

8-711 Range

- Round handle, length 115mm
- Used with DK I/A tips
- I/A tips purchased separately
- Single/Multi thread/Luer Lock options

Ogawa Cannula & Infusion Cannula

8-635-3 18g




- 18 gauge, thin wall, 0.8mm irrigation port, 0.3mm aspiration port
- 30° angled shaft, tip to angle length 7mm
- Coaxial irrigation and aspiration
- Irrigation through main hub and aspiration through side port
- Silicone tubing and luer fitting supplied
- Knurled hub enhances grasp
- Textured tip



Cannula

8-616-1 23g



- 23 gauge, thin wall, 0.7mm diameter x 4mm length
- Straight shaft, 45° bevelled tip
- Grooves in cannula body ensure non-traumatic securing
- Fits through 0.8 to 0.9mm paracentesis opening
- Silicone tubing and luer fitting supplied
- Silicone tubing length 410mm



Infusion Cannula




GEUDER® Stainless Steel-Reusable

Leading manufacturer of high quality surgical instruments for ophthalmology. It has always been the maxim of GEUDER to guarantee the highest degree of precision and reliability.

G22100 D.D. Koch – T. Kohnen Aspiration Handpiece

- 22 gauge 0.7mm, 0.3mm port, textured tip, curved shaft, male connector.




G22101 D.D. Koch – T. Kohnen Irrigation Handpiece

- 22 gauge 0.7mm, 0.7mm irrigation port, lengthened tip, curved shaft.



G15250 Sautter Hydrodissection Cannula

- 27 gauge, 0.4mm, curved shaft.



G32166 Sautter Hydrodissection Cannula

- 27 gauge, 0.45mm.




G32769 H.R. Koch Aspiration Cannula

- 21 gauge, 0.8mm, angled, rough.



G16125 Kratz Capsule Polisher

- Partially etched tip.




G32167 Sautter Hydrodissection Cannula

- 27 gauge, 0.40mm, rotation of 90°.



G33911 Zirm Lasik Irrigating Cannula

- 20 gauge, 0.9mm.






STERIMEDIX® Single-Use-Disposable



Single-use, disposable surgical products meet the demands of modern ophthalmic surgical techniques.

Delivering Global Innovation in Ophthalmic Surgery for More Than 25 Years

STERIMEDIX® strives to continuously improve its service and products to provide innovation and the highest quality, single-use ophthalmic surgical instruments at competitive prices. The entire range of single-use disposable STERIMEDIX® products is delivered in **user-friendly packaging systems** that meet customers' needs and specifications for patient care, inventory control, and operations management.



New! Non-Sterile Bulk Packaging STERIMEDIX® Protector for Kit Packers

Introducing a new, labeled protector for bulk non-sterile products that provides improved quality, flexibility, and reduces waste and cost.

With a specially designed case that protects the cannulas and needle tips, we can ensure that no distortion or deformation will occur to formed or angled products such as cystotomes.

STERIMEDIX® non-sterile protectors are marked with product information and lot traceability nomenclature for effective quality control.

For Inclusion in Kit Packs:

- ✓ Anterior Chamber Maintenance Cannulas
- ✓ Formed Irrigating / Capsularhexis Cystotomes
- ✓ Hydrodissection Cannulas
- ✓ Capsule Polishing Cannulas
- ✓ Refractive / Irrigating Lasik Cannulas
- ✓ I/A (Irrigation / Aspiration) Cannulas



Sterile Packaging Qty 10 / Qty 100 Boxes

STERIMEDIX® user-friendly packaging offers the option of boxes with blister-pack quantities of either 10 or 100 to make it easier to plan and manage stock status for inventory control.

Boxes of 100 offer extra value and help maximize the efficiency of your purchasing dollars.


Keep boxed products on the shelf for backup and unplanned emergencies.



STERIMEDIX® Product Key / Hub Colors

Yellow			30G
Grey			27G
Brown			26G
Orange			25G
Purple			24G
Blue			23G
Black			22G
Green			21G
Yellow			20G
Cream			19G
Pink			18G
White			16G


Anterior Chamber Maintenance Cannulas




Formed Irrigating / Capsularhexis Cystotomes




Hydrodissection Cannulas




Capsule Polishing Cannulas



Refractive / Irrigating Lasik Cannulas



I/A (Irrigation / Aspiration) Cannulas



www.sterimedix.com



STERIMEDIX® Single-Use-Disposable

Single-use, disposable surgical products meet the demands of modern ophthalmic surgical techniques.



STERIMEDIX® Single-Use Product Range

- Refractive Cannulas
- Scleral Markers
- Anesthesia Needles & Cannulas
- Anterior Chamber Cannulas & Maintainers
- Transzonular Injection Cannula
- Washout Cannula
- Viscoelastic Cannulas
- Cystotomes
- Hydrodissection Cannulas
- Viscoexpression Cannulas
- Cortical Cleaving Cannulas
- Nucleus Removal Cannulas
- Cortex Removal Cannulas
- Capsule Polishing Cannulas
- Irrigating/Aspirating Cannulas & Handpieces
- Simcoe Irrigating / Aspirating Cannulas
- Irrigating Microsurgical Handle
- Bimanual Irrigating / Aspirating Handpieces
- Bimanual I/A Systems & Bimanual Cannulas
- Lachrymal Cannulas and DCR Sets
- Vitreoretinal Extension Tube & Self-Retaining Infusion Cannulas
- Trocar Sets
- Vitreoretinal Cannulas
- Vitreoretinal Soft Tips & Brushes
- Backflush and Heavy Liquid Infusion Handles and Combination Sets
- Backflush Handles and Combination Sets
- Charles Flute Handles and Combination Sets
- Non-Sterile Packaging for Kit Packers

STERIMEDIX® strives to continuously improve its service and products to provide innovation and the highest quality, single-use ophthalmic surgical instruments at competitive prices. The entire range of single-use disposable STERIMEDIX® products is delivered in user-friendly packaging systems that meet customers' needs and specifications for patient care, inventory control, and operations management.

One Source®
SERVICE AND SUPPORT

Crestpoint ROUNDS is published monthly by Crestpoint Ophthalmics

Count on us for trusted, personalized service focused on innovation, quality, and comprehensive product selection.

We are the exclusive USA distributors for:

MANI®
DUCKWORTH & KENT®
GEUDER®
STERIMEDIX®

www.sterimedix.com

See You In Chicago

At AAO 2018
October 27-30, 2018

The Crestpoint team will be ready to meet and assist you at AAO 2018 in Chicago! Look for the big blue tower on the main trade show floor!

See you in Chicago

Learn More

ABOUT US

CONTACT US

Mission Statement
Crestpoint Ophthalmics is dedicated to being the foremost, industry leading, single-source provider of the world's highest quality ophthalmic products.

We're your source for 5,000+ hand-held instruments, instrument repair service, repair/replacement program, custom instruments and reusable/disposable instruments. Our top priority is superb customer service to you.

GSA/FSS Contract
V797P-4438B
Ophthalmic Instruments
Surgical Devices
SIN A8a & A8b
Effective
5/1/2011 to 4/30/2021

All information in this publication is subject to change without notice. Trademarks and copyrighted materials are property of their respective owners and are used with permission.

NWBOC
WOSB CERTIFIED: WOMAN OWNED SMALL BUSINESS

All product ranges in this Newsletter are exclusive in the USA.

Duckworth & Kent products are exclusive in the USA and Canada.

Copyright 2018 Crestpoint Ophthalmics. All rights reserved.

Crestpoint Ophthalmics
9100 Watson Road
St. Louis, MO 63126 USA

Phone: (800) 756-7773
Phone: (314) 849-7773
Fax: (314) 849-7776
E-mail: custsvc@crestpointmgt.com
Web: www.crestpointmgt.com