

Crestpoint ROUNDS™

JULY
2018
Volume 2
Issue 6

PUBLISHED FOR ALL OPHTHALMIC AUDIENCES: SURGEONS, YOUNG EYE SURGEONS, NURSES/TECHNICIANS AND PURCHASING TEAMS

IN THIS ISSUE

GEUDER®.....	Cover, 3, 4
Duckworth & Kent®.....	2
Crestpoint Team Spotlight.....	4

Looking forward to seeing you at
AAO 2018
October 27-30 Chicago

MOST POPULAR

GEUDER®
SNYDER-OSHER IOL EXPLANATION SET
SEE ADDITIONAL DETAILS ON PAGE 4

Ophthalmic Surgical Scissors

Ophthalmology is a highly specific surgical discipline that places significant demands on skilled eye surgeons, and their staff. Modern treatments are made possible by precise surgical instruments that meet the exact requirements of the surgical team for quality and reliability. In this mid-year issue of *Crestpoint ROUNDS*, we are focusing on scissors made by the two finest instrument manufacturers in the world: Duckworth & Kent, and GEUDER® AG, and we are proud to say that Crestpoint Ophthalmics is the exclusive North American distributor for both.

Featuring Chromium Plated Stainless-Steel Scissors from GEUDER®


The benefits of GEUDER® chromium plated surgical instruments are:

- ✓ Improved corrosion protection and a professional, high polish finish
- ✓ Ease of cleaning and sterilization helps busy practices keep on schedule
- ✓ Increase in surface hardness and rust prevention from the layer of plating
- ✓ Reduced friction due to the smoother chromium surface
- ✓ Improved liquid flow-off which facilitates the drying process of the instrument



Chromium Plated Stainless-Steel
The sharpest scissors available are made by Geuder from a unique stainless-steel formula. With proper care, these instruments will last a lifetime, and provide exceptional value.

Precision Made In Germany
Geuder has been making surgical instruments for 67 years. Their skilled craftsmen design and build the finest quality, long-lasting, reusable ophthalmic surgical instruments.



Main Instrument Catalog
2018-06



Spring / Summer 2018
Featured Instruments




Geuder Chromium Plating
Certification

MAKE IT A CHROMIUM SUMMER WITH GEUDER®

Now through August 31, 2018 - buy one instrument, select a second at half-price!*





Save with Crestpoint Ophthalmics and GEUDER®
Call (314) 849-7773 or email custsvc@crestpointmgt.com for a quote.
(*equal or lesser value, GEUDER instrument collection)



Duckworth & Kent Titanium-Reusable

Original manufacturer of titanium surgical instruments for ophthalmology and known throughout the world for quality and innovation. D&K Titanium instruments are light weight, non-magnetic and anodized to reduce glare under the microscope.

1-110 VANNAS CAPSULE / CAPSULOTOMY SCISSORS




Vannas Capsule Scissors: Sharp tips, straight blades, tag spring, cut length 5.0mm, flat handle. Titanium, Reusable.

1-111 VANNAS CAPSULE / CAPSULOTOMY SCISSORS




Vannas Capsule Scissors: Sharp tips, curved blades, tag spring, cut length 5.0mm, flat handle. Titanium, Reusable.

1-118 CASTROVIEJO MINIATURE CORNEAL SCISSORS




Castroviejo Miniature Corneal Scissors: Blunt tips, curved blades, cut length 6.0mm, tag spring, flat handle. Titanium, Reusable.

1-218 ANWAR CORNEAL RIGHT SCISSORS




Anwar Corneal Right Scissors: Blunt tips, blades right, polished inner blade, cut length 5.0mm, tag spring, round handle. Titanium, Reusable.

1-219 ANWAR CORNEAL LEFT SCISSORS




Anwar Corneal Left Scissors: Blunt tips, blades left, polished outer blade, cut length 5.0mm, tag spring, round handle. Titanium, Reusable.

1-501 / 1-501B WESTCOTT TYPE STITCH SCISSORS




Westcott Type Stitch Scissors: curved blades, cut length 12.0mm, tag spring, flat handle. Titanium, Reusable.
1-501 semi-sharp tips, curved blades
1-501B blunt tips, curved blades

1-510 WESTCOTT TYPE TENOTOMY SCISSORS



Westcott Type Tenotomy Scissors: Blunt tips, curved blades, cut length 12.0mm, tag spring, flat handle. Titanium, Reusable.

1-705 OSHER IOL CUTTER



Osher IOL Cutter: Angled shafts, cut length 8.0mm, tag spring, flat handle. Micro-serrated cutting edge. Titanium, Reusable.

Geuder
Precision made in Germany

GEUDER[®] Chromium Plated Stainless Steel—Reusable

Leading manufacturer of high quality surgical instruments for ophthalmology. It has always been the maxim of GEUDER to guarantee the highest degree of precision and reliability.


UTILITY RING HANDLE EYE SCISSORS



Curved, pointed-pointed
G-19370 9cm
G-19375 10cm
G-19380 11cm
G-19385 12cm



STRABISMUS SCISSORS



G-19420 straight, long blade, 10.5cm
G-19425 curved, long blade, 10.5cm



STEVENS TENOTOMY SCISSORS



G-19470 straight, short blade, 10.5cm
G-19475 curved, short blade, 10.5cm




ENUCLEATION SCISSORS



G-19520 strong curve, 13.5cm
G-19525 medium curve, 13.5cm
G-19530 light curve, 13.5cm




CASTROVIEJO CORNEAL SCISSORS



G-19690 10.5cm, right
G-19695 10.5cm, left




WESTCOTT SCISSORS



G-19750 blunt-blunt, 12cm
G-19751 pointed-pointed, 12cm



CORNEAL / CORNEAL SECTION SCISSORS



G-19600 light curve, pointed-pointed, 10cm
G-19605 light curve, blunt-blunt, 10cm



VANNAS CAPSULOTOMY SCISSORS



G-19760 extremely delicate, straight, 8cm




VANNAS CAPSULOTOMY SCISSORS




G-19765 extremely delicate, curved, 8cm




HR KOCH CAPSULE SCISSORS




G-32940 45° angled
diameter of shaft 22 gauge / 0.7mm
length of shaft 22mm, handle rotatable



SNYDER-OSHER IOL SCISSORS




For cutting acrylic and silicone foldable IOLs, lower notched blade to create IOL traction, 18 gauge / 1.2mm shaft
G-32996 long model, 4.5mm blades
G-32998 short model, 3.5mm blades
(*see Page 4 for more information)






GEUDER[®] Snyder-Osher IOL Explantation Set



Working successfully through smaller incisions drives this design. The Snyder-Osher scissors have fine serrations that grasp IOL materials during each cut. The forceps have channels on the inside of the rectangular jaws that hold the IOL optic firmly during cutting.






SNYDER-OSHER IOL EXPLANTATION SET

G-32996 Snyder-Osher Scissors

For cutting acrylic and silicone foldable IOLs, lower notched blade to create IOL traction, 18 gauge / 1.2mm shaft, long model, 4.5mm blades (alternate: G-32998 short model, 3.5mm blades)

G-32997 Snyder-Osher Forceps

To hold the IOL while cutting it inside the capsular bag, 2mm long jaws with 0.5mm serrated tips, 22 gauge / 0.7mm lightly curved shaft, overall length 99mm



One Source[®]
SERVICE AND SUPPORT

Crestpoint ROUNDS is published monthly by Crestpoint Ophthalmics

Count on us for trusted, personalized service focused on innovation, quality, and comprehensive product selection.

We are the exclusive USA distributors for:

MANI[®]
DUCKWORTH & KENT[®]
GEUDER[®]
STERIMEDIX[®]

Crestpoint Team Spotlight

Gennie Mason, RN, CRNO

Gennie is one of our expert Customer Resource Reps, and she describes herself as "a nurse helping nurses," because she has worked most of her career as a surgical nurse in the hospital system.

She has a thorough understanding of what hospital nurses need and value when it comes to purchasing, maintaining, and servicing the ophthalmic surgical instruments necessary to keep a busy set of operating rooms running smoothly.

Gennie provides our customers with sales assistance and consulting on anything from a single instrument purchase, to planning and equipping an entire new facility.

See You In Chicago

At AAO 2018
October 27-30, 2018

The Crestpoint team will be ready to meet and assist you at AAO 2018 in Chicago! Look for the big blue tower on the main trade show floor!

See you in Chicago

Learn More

ABOUT US

Mission Statement

Crestpoint Ophthalmics is dedicated to being the foremost, industry leading, single-source provider of the world's highest quality ophthalmic products.

We're your source for 5,000+ hand-held instruments, instrument repair service, repair/replacement program, custom instruments and reusable/disposable instruments. Our top priority is superb customer service to you.

GSA/FSS Contract V797P-4438B
Ophthalmic Instruments Surgical Devices
SIN A8a & A8b
Effective 5/1/2011 to 4/30/2021

All information in this publication is subject to change without notice. Trademarks and copyrighted materials are property of their respective owners and are used with permission.

NWBOC
WOSB CERTIFIED-WOMAN OWNED SMALL BUSINESS

All product ranges in this Newsletter are exclusive in the USA.

Duckworth & Kent products are exclusive in the USA and Canada.

Copyright 2018 Crestpoint Ophthalmics. All rights reserved.

CONTACT US

**Crestpoint Ophthalmics
Crestpoint Management**
9100 Watson Road
St. Louis, MO 63126 USA

Phone: (800) 756-7773
Phone: (314) 849-7773
Fax: (314) 849-7776

E-mail: custsvc@crestpointmgt.com
Web: www.crestpointmgt.com

ZCC483NL (07/2018)