

Crestpoint ROUNDS™

AUGUST
2018
Volume 2
Issue 7

PUBLISHED FOR ALL OPHTHALMIC AUDIENCES: SURGEONS, YOUNG EYE SURGEONS, NURSES/TECHNICIANS AND PURCHASING TEAMS

IN THIS ISSUE


GEUDER®.....	Cover, 2, 4
Duckworth & Kent®.....	3
Crestpoint Team Spotlight.....	4

Newsletter Feature

Nucleus Dividing / Chopping Spatulas / Rotation

There are thousands of types and models of these instruments.

Frequent design changes create an extensive selection available from the best manufacturers.



GEUDER® 2018 Instrument Catalog


We are pleased to announce that GEUDER® has released a comprehensive 2018 instrument catalog with ~1,400 products. This is a well-organized and convenient resource for the ophthalmic surgical supply chain – and it helps us *connect with our customers* to provide sales, service, and support.

✓ Complete reference for GEUDER® surgical instruments

✓ Updated 2018 product selection

✓ Organized with Table of Contents and Indexes

✓ Contact Crestpoint to get your copy



SEE YOU AT **AAO 2018** IN CHICAGO!
EXPERIENCE OUR CUSTOMER SERVICE TEAM **IN PERSON!**



Duckworth & Kent Nucleus Dividers Brochure ZDC427

Featuring detailed information on the procedures and instruments related to:

- > Horizontal, Vertical and Universal Chopping
- > Femto Instruments
- > PreChopping
- > Hooks/Manipulators/Rotators

Includes tip illustrations and examples of use.

[Click to download.](#)







GEUDER[®] Chromium Plated Stainless-Steel. Reusable.


Leading manufacturer of high quality surgical instruments for ophthalmology, made with chromium plated stainless-steel. It has always been the maxim of GEUDER to guarantee the superior level of craftsmanship, precision and reliability.

G16175 DARDENNE IRIS HOOK





DARDENNE Iris hook. Push-pull.
Chromium plated stainless-steel. Reusable.

G16185 DRYSDALE NUCLEUS MANIPULATOR / CAPSULE POLISHER




DRYSDALE Nucleus manipulator and anterior capsule polisher.
Chromium plated stainless-steel. Reusable.

G31507 N AGAHARA NUCLEUS DIVIDING HOOK




N AGAHARA nucleus dividing hook.
Chromium plated stainless-steel. Reusable.

G32505 THOMAS NEUHANN NUCLEUS DIVIDING HOOK





THOMAS NEUHANN nucleus dividing hook.
Chromium plated stainless-steel. Reusable.

G33943 PFAEFFL FEMTO FLAP LIFTER





V-shaped hook with straight front part. Chromium plated stainless-steel. Reusable.

G33950 CHANSUE FEMTO DOUBLE INSTRUMENT



Lenticule spatula, blunt spoon and hook. Chromium plated stainless-steel. Reusable.

G33954 GUELL FEMTO DOUBLE INSTRUMENT




Lenticule spatula, blunt spoon, rounded tip, slightly curved and with hook. Chromium plated stainless-steel. Reusable.

G38607 DSAEK DESCemet INCISION HOOK



Diameter of hook: 0.2mm
Total length: 103mm
Chromium plated stainless-steel. Reusable.




Duckworth & Kent Titanium. Reusable.

Original manufacturer of titanium surgical instruments for ophthalmology and known throughout the world for quality and innovation. D&K Titanium instruments are light weight, non-magnetic and anodized to reduce glare under the microscope.



Choppers

6-085-1 Barrett Duo Nucleus Rotator / Manipulator / Splitter



Nucleus splitter

- Straight, sharp inner / outer sides of tip, tip length 1.25mm
- Cutting edges 60° to axis
- 45° angled shaft, tip to angle length 14mm



Rotator / manipulator


- 0.65mm mushroom tip
- 45° angled shaft, tip to angle length 10mm
- Round handle, length 130mm

Manipulators


6-400 Kuglen Iris Hook and IOL Manipulator




- 0.76mm x 0.75mm x 0.25mm rotator tip
- Push, pull style
- Angled shaft, tip to angle length 9mm
- Round handle, length 122mm



6-496 Drysdale Rotator




- 0.5mm tip with 1.3mm paddle
- 45° angled shaft, tip to angle length 10mm
- Round handle, length 118mm




Hooks

6-250 DK Sinskey Hook




- 0.18mm diameter tip pointing down
- Angled shaft, tip to angle length 10mm
- Round handle, length 119mm




Probes

6-656 DK Lacrimal Probe




- High quality surface finish over first 5mm of tips
- Overall length 130mm
- Tip sizes: 0000 and 000




Spatulas

6-870 Buratto LASIK Oval Spatula




- 1.3mm curved blade
- 35° angled curved shaft, tip to angle length 10.5mm
- Round handle, length 122mm



Designed to raise flap during LASIK.

CORE VALUES Collection

Duckworth & Kent is proud to present the Core Values Collection a range of reusable instruments offering our renowned precision and quality at an exceptional price.



See leaflet for full range of instruments

Visit www.duckworth-and-kent.com to view our 'Core Values Collection' and much more.



GEUDER® Summer Sale Extended! September 30, 2018!

Summer isn't over yet! **By popular demand**, we have extended our GEUDER Summer Sale* now to September 30, 2018, buy one GEUDER instrument, select a second GEUDER instrument at half-price! (equal or lesser value, GEUDER instrument collection)

MAKE IT A CHROMIUM SUMMER WITH GEUDER®
 Extended!! Now through September 30, 2018 — buy one instrument, select a second at half-price*

GEUDER® Builds the Sharpest Scissors
 Known worldwide for the sharpest scissors. Made in Germany by skilled craftsmen.

Chromium-Plated Stainless-Steel
 Adds corrosion resistance and hardness to the surface, making the instrument easier to clean, sterilize, and maintain.

Precision Engineering and Machining
 Geuder's team of 40 surgical, microsurgical, and precision instrument engineers are among the best in their field.

Save With Crestpoint Ophthalmics and GEUDER®
 It's the perfect season to acquire Geuder chromium-plated instruments to improve surgical operations and patient outcomes. Call (314) 849-7773 or email custsvc@crestpointmgt.com for a quote. (*equal or lesser value, GEUDER instrument collection)

*May not be combined with other offers. Contact Crestpoint to receive your copy of the 2018 GEUDER Instrument Catalog.



One Source®
 SERVICE AND SUPPORT

Crestpoint ROUNDS is published monthly by Crestpoint Ophthalmics

Count on us for trusted, personalized service focused on innovation, quality, and comprehensive product selection.

We are the exclusive USA distributors for:
 DUCKWORTH & KENT®
 GEUDER®
 MANI®

We are the USA distributors for:
 STERIMEDIX® | Q INSTRUMENTS®

Crestpoint Team Spotlight

Gayle Sindel

Gayle manages our accounts payable and accounts receivable department and she's a key member of our operations team, helping to make your billing and payment experience with Crestpoint smooth and seamless.

A seven year veteran of our firm, Gayle says she "has the best job in the world" because she gets to interact with our customers and vendors on a daily basis. At home, Gayle is an avid reader, and she's looking forward to welcoming her first grandchild in April 2019.

We always strive to make your experience the best it can be, and Gayle is one of our helpful ambassadors. We can't thank our customers enough!

See You In Chicago At AAO 2018
 October 27-30, 2018

Crestpoint / D&K / GEUDER
 Booths #4200 & #3900
 MANI #4203
 STERIMEDIX #4303
 AMCON #3304

See you in Chicago [Learn More](#)

ABOUT US

Mission Statement
 Crestpoint Ophthalmics is dedicated to being the foremost, industry leading, single-source provider of the world's highest quality ophthalmic products.

We're your source for 5,000+ hand-held instruments, instrument repair service, repair/replacement program, custom instruments and reusable/disposable instruments. Our top priority is superb customer service to you.

Crestpoint OPTHALMICS

GSA/FSS Contract V797P-4438B
 Ophthalmic Instruments Surgical Devices
 SIN A8a & A8b Effective 5/1/2011 to 4/30/2021

All information in this publication is subject to change without notice. Trademarks and copyrighted materials are property of their respective owners and are used with permission.

NWBOC
 WWSB CERTIFIED-WOMAN OWNED SMALL BUSINESS

All product ranges in this Newsletter are exclusive in the USA.

Duckworth & Kent products are exclusive in the USA and Canada.

Copyright 2018 Crestpoint Ophthalmics. All rights reserved.

CONTACT US

**Crestpoint Ophthalmics
 Crestpoint Management**
 9100 Watson Road
 St. Louis, MO 63126 USA

Phone: (800) 756-7773
 Phone: (314) 849-7773
 Fax: (314) 849-7776
 E-mail: custsvc@crestpointmgt.com
 Web: www.crestpointmgt.com

ZCC483NL (08/2018)