

# Crestpoint ROUNDS™

SEPTEMBER  
2017  
Volume 1  
Issue 1

PUBLISHED FOR ALL OPHTHALMIC AUDIENCES: SURGEONS, YOUNG EYE SURGEONS, NURSES/TECHNICIANS AND PURCHASING TEAMS

## IN THIS ISSUE

Six New Catalogs Available....	1
Specula and Side Port Incisions:	
<i>Duckworth &amp; Kent</i> .....	2
<i>MANI</i> .....	2
<i>Geuder</i> .....	3
<i>Sterimedix</i> .....	3
<i>Mr.Phaco/Mr.Femto</i> .....	4
Our Office on Route 66.....	4
AAO New Orleans.....	4
About Us.....	4
Contact Us.....	4

## Specula and Side Port Incisions

Crestpoint is proud to introduce *Crestpoint Rounds*. Published monthly, each Newsletter will highlight featured products. Please contact us to discuss additional product styles available from over 5,000+ items in our full product range.

In this issue, Specula and Side Port incisions.

See you at AAO 2017 in New Orleans, Booths #816, #919, #921 and #925.

We cannot thank our customers enough!



**J. Michelle Glossip**  
Chief Executive



**Henry M. Mitchell**  
President




**Dot C. Smith**  
Customer Resource  
Director




**Michael Dill**  
Customer Resource  
Manager

## Six New Catalogs Available

We strive to provide our customers with **innovative products** and the **finest surgical instruments available**. Extensive product lines for every step of your procedures. E-mail or call to request your copies. Click below for more information.






# Duckworth & Kent

Original manufacturer of titanium surgical instruments for ophthalmology.

### Kratz Barraquer Speculum

9-555 Straight  
9-556 Temporal  
9-556F Temporal (heavy)  
9-557 Nasal  
**Features**  
Single piece construction, open blades, 14.5mm wide


Titanium, Reusable



### Lieberman Adjustable Speculum

9-588-3 Temporal  
**Features**  
Temporal, open blades, adjustable with thumb screw, 15mm wide


Titanium, Reusable



### Fine Thornton Fixation Ring

9-510 Flat Handle  
9-510R Round Handle  
**Features**  
3/4, 13.0mm, 9 point fixation, pivot swivel, flat handle


Titanium, Reusable



### Shepard Reversible Speculum


9-599 Solid Blades  
9-599-1 Closed Blades  
**Features**  
14mm, reversible nasal or temporal

Titanium, Reusable



### Single Piece


All D&K Specula are made from a single piece of titanium. There are no joints which create weak points in the construction.



### Adjustable

Blades can open to 30mm (without closing resistance)

Main body is made from a single piece of titanium with added strength in the arms to give a rigid strong construction, which overcomes any attempt by patient to squeeze on the blades.



Blades are reduced in thickness to 0.5mm. This gives added comfort to the patient, yet the strength is not compromised.

Smooth, light action adjustable thumb screw mechanism with large easy to grip knurled thimble.



# MANI®

Specializes in single use Knives, Sutures, Trocars for all ophthalmic procedures. **New! Marking Straight Knife.**

### Straight Knife

MST15  
**Features**  
Main Application: side port, 6pcs/box


Disposable




### Marking Straight Knife

MST15MK  
**Features**  
Main Application: side port, 6pcs/box. By placing the lines serving as guides of the side port width in 1.0mm and 1.2mm apart, the operability has been further improved.

Disposable



MANI's fine-processing technique has created the crosslapped edge, which provides excellent sharpness and consistent quality. Color-coded handles allow for easy identification. The color is different depending on the tip shape.



MANI slit knife - front side with Non-Glare Coating

The crosslapped edge made of MANI's original hard fiber stainless steel produces the ultimate sharpness. In order to maintain the sharpness, a patented Non-Glare treatment is processed on knives except its edge and back side.






# Geuder®

Leading manufacturer of high quality surgical instruments for ophthalmology. It has always been the maxim of Geuder to guarantee the highest degree of precision and reliability.

### Schott Adjustable Speculum

**G15921 Adult**  
**G15926 Pediatric**  
**G15928 Pediatric/Aspiration**  
**Features**  
 Solid blades, adjustable


Stainless Steel, Reusable



### Leiberman-Tennant Adjustable Speculum

**G15960 Temporal**  
**G15961 Temporal/Aspiration**  
**Features**  
 Fenestrated, open blades for temporal placement, adjustable


Stainless Steel, Reusable



### Kershner Adjustable Speculum

**G15965 Nasal & Temporal**  
**Features**  
 Fenestrated, closed blades, for nasal and temporal placement, adjustable


Stainless Steel, Reusable



### Solomon Adjustable Speculum

**G33971**  
**Features**  
 Small, open valves with middle bar, double knurled screw for fine tuning, adjustable, 14.5mm wide


Stainless Steel, Reusable



### Dick Adjustable Speculum

**G33965**  
**Features**  
 Open valves with middle bar, adjustable, 18mm wide

Stainless Steel, Reusable



### Cook Adjustable Speculum

**G15995**  
**Features**  
 Medium, solid blades, adjustable

Stainless Steel, Reusable



# Sterimedix®

Complete range of single use devices for all ophthalmic surgery procedures.

<b>M2639</b>	25G x 1"	0.5 x 25mm	Anterior Chamber Cannula. Angled 35°, 7.5mm from tip (Rycroft)
<b>M2270</b>	30G x 7/8"	0.3 x 22mm	Anterior Chamber Cannula. Angled 35°, 4mm from tip (Rycroft)
<b>M2270A</b>	30G x 7/8"	0.3 x 22mm	Anterior Chamber Cannula. Angled 35°, 8mm from tip (Rycroft)
<b>M2273</b>	30G x 7/8"	0.3 x 22mm	Anterior Chamber Cannula. Angled 35°, 5.5mm from tip (Rycroft)
<b>M2273A</b>	27G x 7/8"	0.4 x 22mm	Anterior Chamber Cannula. Angled 35°, 4.5mm from tip (Rycroft)
<b>M2279</b>	27G x 7/8"	0.4 x 22mm	Anterior Chamber Cannula. Angled 35°, 11mm from tip (Rycroft)
<b>M2275</b>	25G x 7/8"	0.5 x 22mm	Anterior Chamber Cannula. Angled 35°, 5mm from tip (Rycroft)
<b>M2273G</b>	23G x 7/8"	0.64 x 22mm	Anterior Chamber Cannula. Angled 35°, 5.5mm from tip (Rycroft)
<b>M2642</b>	21G x 7/8"	0.81 x 22mm	Anterior Chamber Cannula. Angled 35°, 7.5mm from tip (Rycroft)
<b>M2273B</b>	20G x 7/8"	0.9 x 22mm	Anterior Chamber Cannula. Angled 35°, 7.5mm from tip (Rycroft)





# Mr.Phaco® / Mr.Femto®

One Source service and support for your ophthalmic needs. New innovative instruments, repair/replacement program, custom instruments. Excellent reputation for service and quality.

### Thornton Davidson Fixation Ring

Stainless Steel, Reusable

#### PH2115 Flat Handle

**Features**  
3/4, pivot swivel, flat handle, bead blasted double sided finish for good grip, no fixation points, gentle on conjunctiva



### Modi Speculum

Stainless Steel, Reusable

#### PH9004

**Features**  
Temporal, open blades, adjustable, Silastic tubing on blades, gentle on conjunctiva



## YOUNG EYE SURGEONS PROGRAM

*Crestpoint offers assistance to the new surgeon during the transition period from training to practice. We can help with choosing product, selecting tip and handle styles, and by offering special pricing on purchases. Please contact Crestpoint for more information.*

### Route 66

Our offices are located on Watson Road in Crestwood, Missouri — and our address has a lot of history behind it. Watson is the modern day Route 66, a highway that changed the landscape of America. Drive two miles northeast of our offices on Watson and you'll find the Wayside Motel, which was built in 1940, and still has its original carports. Keep going another two miles and you'll end up at Ted Drewes Frozen Custard, where you can treat yourself to the best dessert in the world (no, we're not biased). "From Chicago to L.A.," there's so much to explore on Route 66, and we're proud to be a part of this slice of American history.

### AAO 2016 Chicago



### AAO 2017

**SEE YOU IN NEW ORLEANS**  
November 11-14, 2017  
Booths #816, #919, #921, #925

Crestpoint is gearing up for AAO 2017 in New Orleans! We're making some changes to our booth and can't wait to see you at the show!

Look for our brilliant blue tower, and our extensive product display, featuring the latest. Watch for more information in our next month's *Crestpoint Rounds*.

One Source®

## ABOUT US

### Mission Statement

Crestpoint Ophthalmics is dedicated to being the foremost, industry leading, single-source provider of the world's highest quality ophthalmic products.

We're your source for 5,000+ hand-held instruments, instrument repair service, repair/replacement program, custom instruments and reusable/disposable instruments. Our top priority is superb customer service to you.



All information in this publication is subject to change without notice.

Trademarks and copyrighted materials are property of their respective owners and are used with permission.

**NWBOC**  
WOSB CERTIFIED-WOMAN OWNED SMALL BUSINESS

All product ranges in this Newsletter are exclusive in the USA.

## CONTACT US

### Crestpoint Management

9100 Watson Road  
St. Louis, MO 63126 USA

Phone: (800) 756-7773

Phone: (314) 849-7773

Fax: (314) 849-7776

E-mail: [custsvc@crestpointmgt.com](mailto:custsvc@crestpointmgt.com)

Web: [www.crestpointmgt.com](http://www.crestpointmgt.com)

ZCC483NL (09/2017)