

Crestpoint ROUNDS™

SEPTEMBER
2018
Volume 2
Issue 8

PUBLISHED FOR ALL OPHTHALMIC AUDIENCES: SURGEONS, YOUNG EYE SURGEONS, NURSES/TECHNICIANS AND PURCHASING TEAMS

IN THIS ISSUE

Duckworth & Kent®	Cover, 2, 4
GEUDER®	3
Crestpoint Team Spotlight	4

MANUFACTURER FEATURE

GUIDANCE ON THE WORLD'S BEST OPHTHALMIC SURGICAL INSTRUMENTS FOR YOUR PRACTICE

Duckworth & Kent is the undisputed world leader in the design, engineering, and production of titanium alloy precision ophthalmic surgical instruments.

With increasing technological challenges and demanding customer expectations, the value of innovative design and advanced manufacturing is stronger now than ever before.

By exploring advanced concepts and embracing new technology, Duckworth & Kent remains at the forefront of ophthalmic surgical instrument design, engineering, and quality.




Precision Engineered Single Piece, and Adjustable Specula

Featuring long-lasting reusable titanium instruments from Duckworth & Kent®

Why is titanium alloy important?

- ✓ Lightweight and extremely durable
- ✓ Corrosion resistant and no oxidization
- ✓ Non-magnetic, non-reflective
- ✓ Withstands repeat sterilization

Non-Reflective Finish
Anodized surface.



Blade Thickness Reduced to 0.5mm
Thinner blades provide additional comfort to the patient while overall strength and shape are maintained.

Smooth, Light Action Adjustable Thumbscrew
Large knurled handle is easy to grip and turn.


D&K Item 9-579-2
Mackkool® Open Blade Speculum

Main Body Single Piece of Titanium
With added strength in the arms to produce a rigid, strong, durable mechanism that holds position.

SEE YOU AT **AAO 2018** IN CHICAGO!
OCTOBER 27-30, BOOTHS #3900 AND #4200



Duckworth & Kent
2018 Catalog



Duckworth & Kent Titanium. Reusable.


Original manufacturer of titanium surgical instruments for ophthalmology and known throughout the world for quality and innovation. D&K Titanium instruments are light weight, non-magnetic and anodized to reduce glare under the microscope.

Specula—Single Piece (many more styles available, contact Crestpoint customer service)

9-551 BARRAQUER ADULT SPECULUM, TEMPORAL

Essential CATARACT SET

- 14.5mm closed blades
- Angled to rest temporally
- Single piece construction
- Titanium, reusable




9-556 KRATZ BARRAQUER SPECULUM, TEMPORAL

Essential CATARACT SET


CORE VALUES Collection

- 14.5mm open blades
- Angled to rest temporally
- Single piece construction
- Titanium, reusable




9-552 BARRAQUER ADULT SPECULUM, NASAL

- 14.5mm closed blades
- Angled to rest nasally
- Single piece construction
- Titanium, reusable



9-572 BARRAQUER PAEDIATRIC SPECULUM, TEMPORAL

- 6mm closed blades
- Straight to rest temporally
- Single piece construction
- Titanium, reusable



Specula—Adjustable (many more styles available, contact Crestpoint customer service)

9-588-2 HORN ADJUSTABLE FEMTOSECOND LASER SPECULUM, TEMPORAL


- 15.5mm open blades
- Curved to rest temporally
- Adjustable thumb screw
- Suitable for femtosecond laser, and can be used for cataract surgery and LASIK
- Blades allow clearance for docking devices
- Rigid design
- Titanium, reusable



9-588-3 LIEBERMAN ADJUSTABLE SPECULUM, TEMPORAL

Essential CATARACT SET

- 15mm open blades
- Angled to rest temporally
- Adjustable thumb screw
- Rigid design
- Titanium, reusable



9-590 BARRETT ADJUSTABLE SPECULUM, TEMPORAL

Essential CATARACT SET

- 14.5mm open blades
- Angled to rest temporally
- Adjustable thumb screw
- Rigid design
- Titanium, reusable



9-598-2 KOCH CIONI OPEN BLADE ADJUSTABLE SPECULUM, NASAL

CORE VALUES Collection

- 13.5mm open blades
- Angled to rest nasally
- Adjustable thumb screw
- Rigid design
- Titanium, reusable



9-595 MACKOOL® ADJUSTABLE SPECULUM, TEMPORAL

- 14.5mm solid blades
- Angled to rest temporally
- Adjustable thumb screw
- Rigid design
- Not compressible
- Critical for LASIK
- Titanium, reusable



9-599-1 SHEPARD REVERSIBLE SPECULUM, TEMPORAL OR NASAL

- 14.0mm closed blades
- Angled to rest temporally or nasally
- Adjustable thumb screw
- Rigid design
- Titanium, reusable



Geuder
Precision made in Germany


GEUDER[®] Chromium Plated Stainless-Steel. Reusable.

Leading manufacturer of high quality surgical instruments for ophthalmology, made with chromium plated stainless-steel. It has always been the maxim of GEUDER to guarantee the superior level of craftsmanship, precision and reliability.

Specula—Wire and Adjustable (many more styles available, contact Crestpoint customer service)


G-32133 KRATZ-BARRAQUER EYE SPECULUM, FENESTRATED

- 15.0mm round wing blades
- Length 40.0mm
- To lift and retract the lids away from the globe, eliminating external pressure
- Chromium plated stainless-steel
- Reusable



G-17004 H.R. KOCH TEMPORAL APPROACH EYE SPECULUM

- 15.0mm round blades
- Length 40.0mm
- For nasal placement
- Fenestrated
- Chromium plated stainless-steel
- Reusable



G-32119 KRATZ-BARRAQUER ADJUSTABLE EYE SPECULUM

- 15.0mm round wing blades
- Length 40.0mm
- Adjustable
- Chromium plated stainless-steel
- Reusable



G-17030 SAUER ADJUSTABLE EYE SPECULUM

- 12mm solid blades
- Length 30mm
- Adjustable
- Chromium plated stainless-steel
- Reusable



Specula—Adjustable / Locking (many more styles available, contact Crestpoint customer service)

G-15946 MELLINGER-BEARD ADJUSTABLE / SELF RETAINING SPECULUM

- Large
- 15.5mm solid blades
- Length 65.0mm
- Adjustable
- Self retaining
- Chromium plated stainless-steel
- Reusable



G-16020 CASTROVIEJO ADJUSTABLE / LOCKING SPECULUM

- 14.0mm blades-small
- Length 96.0mm
- Fenestrated
- Adjustable
- Locking
- Chromium plated stainless-steel
- Reusable




G-15990 COOK ADJUSTABLE / LOCKING SPECULUM

- Large
- 16.0mm solid blades
- Length 50.0mm
- Adjustable
- Locking screw
- Chromium plated stainless-steel
- Reusable




G-15960 LIEBERMAN-TENNANT ADJUSTABLE / LOCKING SPECULUM

- 15.0mm open blades
- Length 80.0mm
- Adjustable
- Fenestrated for temporal placement
- Locking screw
- Chromium plated stainless-steel
- Reusable



G-15966 KERSHNER ADJUSTABLE / REVERSIBLE LOCKING SPECULUM


- Adjustable / reversible
- 15.0mm solid blades
- Length 75.0mm
- For nasal and temporal placement
- Locking screw
- Chromium plated stainless-steel
- Reusable



G-33971 SOLOMON FEMTOCATARACT EYE SPECULUM

- 14.5mm flattened fenestrated blades with middle bar
- Length 90.0mm
- Adjustable
- Double knurled screw for fast adjustment
- Chromium plated stainless-steel
- Reusable





Duckworth & Kent Titanium. Reusable.

Cost effective and eco-friendly. When you invest in a Duckworth & Kent instrument, you are purchasing a product that will serve you and your practice impeccably for many years and countless procedures. D&K's titanium alloy surgical instruments are lightweight, strong and durable, non-magnetic, and corrosion resistant.

Essential CATARACT SET



The decision to choose D&K reusable instruments over single-use disposables is a simple one. When you invest in a set of D&K's high quality, word-renowned titanium alloy instruments, your practice will:

- ✓ Increase efficiency
- ✓ Save money long term
- ✓ Reduce waste and lower your environmental impact

An increase in quality and efficiency means a smoother running organization. But, the most important result is **improved patient care**. Contact Crestpoint Ophthalmics to see how much you can save by investing in D&K's high quality titanium alloy surgical instruments.



One Source
SERVICE AND SUPPORT

Crestpoint ROUNDS is published monthly by Crestpoint Ophthalmics

Count on us for trusted, personalized service focused on innovation, quality, and comprehensive product selection.

We are the exclusive USA distributors for:
DUCKWORTH & KENT®
GEUDER®
MANI®

We are the USA distributors for:
STERIMEDIX® | Q INSTRUMENTS®

Crestpoint Team Spotlight

Dot Smith

Dot is our Customer Resource Director, and she's a key part of our management team, helping to ensure that your customer service experience with Crestpoint is always handled in a positive, efficient, and friendly manner.

Dot was the first employee to join our team when we opened in 1994, and she's now in her twenty-fourth year. She says the best parts of working in customer service are the relationships: "I love working with people, and I never really meet a stranger!"

We always strive to make your experience the best it can be, and Dot is one of our helpful ambassadors. We can't thank our customers enough!

See You In Chicago At AAO 2018
October 27-30, 2018

Crestpoint / D&K / GEUDER
Booths #3900 and #4200
MANI #4203
STERIMEDIX #4303
AMCON #3304

See you in Chicago [Learn More](#)

ABOUT US

CONTACT US

Mission Statement
Crestpoint Ophthalmics is dedicated to being the foremost, industry leading, single-source provider of the world's highest quality ophthalmic products.

We're your source for 5,000+ hand-held instruments, instrument repair service, repair/replacement program, custom instruments and reusable/disposable instruments. Our top priority is superb customer service to you.

Crestpoint OPTHALMICS

GSA/FSS Contract V797P-4438B
Ophthalmic Instruments Surgical Devices
SIN A8a & A8b Effective 5/1/2011 to 4/30/2021

All information in this publication is subject to change without notice. Trademarks and copyrighted materials are property of their respective owners and are used with permission.

NWBOC
WOSB CERTIFIED - WOMAN OWNED SMALL BUSINESS

All product ranges in this Newsletter are exclusive in the USA.

Duckworth & Kent products are exclusive in the USA and Canada.

Copyright 2018 Crestpoint Ophthalmics. All rights reserved.

**Crestpoint Ophthalmics
Crestpoint Management**
9100 Watson Road
St. Louis, MO 63126 USA

Phone: (800) 756-7773
Phone: (314) 849-7773
Fax: (314) 849-7776
E-mail: custsvc@crestpointmgt.com
Web: www.crestpointmgt.com