


Crestpoint ROUNDS™

OCTOBER
2018
Volume 2
Issue 9

PUBLISHED FOR ALL OPHTHALMIC AUDIENCES: SURGEONS, YOUNG EYE SURGEONS, NURSES/TECHNICIANS AND PURCHASING TEAMS

Special Edition—See You In Chicago!

WELCOME TO CRESTPOINT—D&K—GEUDER—MANI—MR.PHACO/MR.FEMTO—STERIMEDIX—AMCON



DUCKWORTH & KENT
E Range Instruments
Page 2

Original manufacturer of titanium instruments for ophthalmic surgery.




GEUDER
New 2018 Catalog
Page 3

High quality stainless steel surgical instruments for ophthalmic surgery.



MANI
Short Head Knives
Page 2

MANI's fine processing technique creates the crisslapped edge (patented).



MR.PHACO
Innovative Products
Page 4

Innovative and custom instruments, repair/replacement program.



STERIMEDIX
New 2018 Catalog
Page 3

Single-use disposable products deliver innovation in ophthalmic surgery.



Research & Development



Mr.Phaco PH6018 Ogawa Stamer Hook
design.

Crestpoint Innovation
Bring your ideas to us, basic or advanced.

We will treat your idea as a priority and innovate to create a new design, or improve the features and characteristics of an existing design.

We're an industry leader with a sterling track record. Combining years of industry experience with significant clinical and medical expertise, our team is uniquely qualified to consult and assist with surgical instrument design, engineering, and production. We have helped design over one thousand ophthalmic surgical devices.


We work directly with surgeons and the world's best manufacturers to research, develop, and produce new and innovative ophthalmic surgical instruments that deliver effective solutions for improved patient care.

Crestpoint at AAO 2018


Experience Crestpoint's largest show of products in the world at #3900, 4200, 4203, 4303, 3304. Book a private appointment to preview any special items that you would like to see. (See floor map Page 4)


The owner, exclusive supply partners, new representatives, and our specialized customer resources team will be on hand. Dedicated to the ophthalmic community, with **best-of-the-best products and service.** Always count on Crestpoint for trusted, personalized attention.

OneSource® for 5,000+ hand-held surgical instruments, sales, service, repair/replacement programs, research & development, design, engineering, manufacturing, and expert technical support. **Our top priority is superb customer service to you.**



AAO 2018 Floor Map (large view Page 4)






Duckworth & Kent Titanium. Reusable.

Original manufacturer of titanium surgical instruments for ophthalmology and known throughout the world for quality and innovation. D&K titanium instruments are lightweight, non-magnetic, and anodized to reduce glare under the microscope.

E-RANGE Instruments




The E-Range represents the application of cutting edge technology fused with over 50 years of engineering knowledge. These instruments offer surgeons an elegant, well balanced and ergonomically designed handle while still providing our high quality, consistent and durable functionality.




2-716GE-1 Inamura Flat Handle Capsulorhexis Forceps

- Sharp, pointed, serrated tips
- Cross action tips, 1.5mm wide at pivot box
- Curved shaft, flat handle
- Tip to pivot 10.5mm
- Tips angled 45° from shaft
- Overall length 111.0mm
- Marks on shaft 2.5mm and 5.0mm denote desired size of capsulorhexis
- Titanium, reusable



2-505ER8 Curved Tying Round Handle Forceps

- 6.0mm tying platforms
- Curved shafts
- Round handle
- Overall length 115.0mm
- Titanium, reusable



E-Range Features

- Ergonomically Designed**
- Expertly Crafted**
- Engineered to Perfection**
- Exceptional Price**

The E-Range is single piece in design (excluding cross action forceps), removing the variability and physical weaknesses associated with instruments produced from component parts. This balanced and creative style of instrument is produced in the most efficient manner, enabling Duckworth & Kent to offer the surgeon an incredible range of instruments at an exceptional price.



MANI® Ophthalmic Knives

MANI's bevel shape (patented), crosslapped edge, and high-strength stainless steel construction produces the "Ultimate Sharpness" for a true, stable incision with higher accuracy and lower penetration resistance. Color-coded knife handles allow for easy identification. The color is different depending on the tip shape.

Short-Head Slit Knife




Main Application: Corneal Incision, Scleral Incision, 6pcs/box


- | | |
|-----------------------------|---|
| MSL20SH – Slit Angled 2.0mm | MSL26SH – Slit Angled 2.65mm |
| MSL21SH – Slit Angled 2.1mm | MSL27SH – Slit Angled 2.75mm |
| MSL22SH – Slit Angled 2.2mm | MSL28SH – Slit Angled 2.8mm |
| MSL23SH – Slit Angled 2.3mm | MSL29SH – Slit Angled 2.9mm |
| MSL24SH – Slit Angled 2.4mm | MSL30SH – Slit Angled 3.0mm |
| MSL25SH – Slit Angled 2.5mm | (photo below—individual safety packaging) |



Example Use of the Marking Lines




Extreme sharpness, simple operation, high visibility. The position of the knife angle is moved 2.0mm towards the tip in order to make the blade smaller while maintaining the same "Ultimate Sharpness" as our conventional models. This design is highly effective for patients with sunken/deep-set eyes or narrow lid fissures. Lower risk of contact with specula or drapes. Laser marking lines at 1.5mm and 2.0mm from the blade tip enable simple operation for creating a tunnel. MSL20~24SH=(8.0mm→6.0mm)/MSL25~30SH=(9.0mm→7.0mm)



Click for video



Click for video




Click for video



GEUDER® Chromium Plated Stainless Steel. Reusable.

Leading manufacturer of high quality surgical instruments for ophthalmology, made with chromium plated stainless steel. It has always been the maxim of GEUDER to guarantee the superior level of craftsmanship, precision and reliability.

New! Geuder 2018 Catalog and Scissors Mini-Brochure



The Geuder 2018 catalog is a comprehensive 312 page reference containing 1,400 instruments, specifications, tip illustrations, and instrument sets. Search features include a detailed table of contents, numerical index, keyword index, and name index, available in both PDF and print versions.

The Crestpoint Geuder Scissors Mini-Brochure is a quick reference guide to the most popular items.

G-22100 D.D. Koch – T. Kohnen Aspiration Handpiece

- 22 gauge / 0.7mm
- 0.3mm aspiration port
- Textured tip
- Chromium plated stainless-steel, reusable
- G-22102 23 gauge / 0.6mm



G-22101 D.D. Koch – T. Kohnen Irrigation Handpiece

- 22 gauge / 0.7mm
- Front opening
- Chromium plated stainless-steel, reusable



G-32996 Snyder – Osher IOL Scissors



- For cutting acrylic and silicone foldable IOLs, lower notched blade to create IOL traction
- 18 gauge / 1.2mm shaft, overall length 150.0mm, chromium plated stainless steel, reusable
- Long model, 4.5mm blades
- G-32998 short model, 3.5mm blades
- Use with G-32997 Snyder-Osher Forceps



STERIMEDIX® Single-Use Disposable

Safety. Quality. Flexibility. Cost efficiency. STERIMEDIX single-use, disposable surgical products meet the demands of modern ophthalmic surgical techniques.

Delivering Global Innovation in Ophthalmic Surgery for More Than 25 Years

STERIMEDIX strives to continuously improve its service and products to provide innovation and the highest quality, single-use ophthalmic surgical instruments at competitive prices. The entire range of single-use disposable STERIMEDIX products is delivered in user-friendly packaging systems that meet customers' needs and specifications for patient care, inventory control, and operations management.

Non-Sterile Bulk Packaging for Kit Packers Sterile Packaging Qty 10 / Qty 100 Boxes



Labeled protector for bulk non-sterile products that provides improved quality, flexibility, and reduces waste and cost. A specially designed case protects the cannulas and needle tips, to ensure no distortion or deformation occurs to formed or angled products such as cystotomes.

Each protector is marked with product details and nomenclature for traceability and quality control.

STERIMEDIX user-friendly packaging offers the option of boxes with blister pack quantities of either 10 or 100 to make it easier to plan and manage stock status for inventory control.



New Products for 2018

Contact us for the new STERIMEDIX 2018 Catalog with full range of products and details, specifications and index.



STERIMEDIX®

Continuous development from Sterimedix

to meet the demands of modern ophthalmic surgery techniques

www.sterimedix.com

STERIMEDIX®

The leading European-based manufacturer of single-use ophthalmic cannulas and handpieces

www.sterimedix.com



Mr.Phaco® / Mr.Femto®

OneSource® service and support for your ophthalmic needs. New innovative instruments, R&D, design, engineering, repair/replacement program, custom instruments. Excellent reputation for service and quality.

MMP265 Osher Sweep



- Sweep shaped similar to traditional muscle hook, but shorter with a different bend to help push escaped aqueous from conjunctiva.
- Stainless steel, reusable

CLICK FOR VIDEO



PH9004 Modi Femto Speculum



- Temporal
- Open blades
- Adjustable thumb screw
- 15.5mm wide
- Non-glare finish
- Stainless steel, reusable

This speculum is ergonomically correct, allowing for backward pressure on globe to effect a good dock.




See You in Chicago at AAO 2018!

Best-of-the-best reusable and disposable products—**Duckworth & Kent, Geuder, MANI, Sterimedix, AMCON, Q** instruments, **Mr.Phaco/Mr.Femto**.

OneSource® for innovation, service, technical support, and repair/replacement/refurbish programs.

Please contact our Customer Resource Team at (314) 849-7773.



Contact Us

Crestpoint Ophthalmics is dedicated to being the foremost, industry leading, single-source provider of the world's highest quality ophthalmic products.

Our top priority is superb customer service to you.



Crestpoint Management
Crestpoint Ophthalmics
 9100 Watson Road
 St. Louis, MO 63126 USA

GSA/FSS Contract
 V797P-4438B
 Ophthalmic Instruments
 Surgical Devices
 SIN A8a & A8b
 Effective
 5/1/2011 to 4/30/2021

Published by Crestpoint Management. Manufacturer trademarks, part data, and copyrighted materials are property of their respective owners and are used with permission. D&K, Geuder, MANI product ranges in this publication are exclusive in the USA / NA.

- Phone: (800) 756-7773
- Phone: (314) 849-7773
- Fax: (314) 849-7776
- E-mail: custsvc@crestpointmgt.com
- Web: www.crestpointmgt.com

Copyright 2018 Crestpoint Management / Crestpoint Ophthalmics
 All rights reserved. All information in this publication is subject to change without notice.

ZCC483NL (10/2018)

